

CREATINA MICROSULES
Amino acids supplement
For the development and strengthen of muscle
Oral powder



Product basic composition

Creatine Monohydrate, excipients q.s.

Guarantee levels per kg of product

Creatine Monohydrate 800 g
Excipients q.s.ad..... 1000 g

Each 100 ml contains

Creatine Monohydrate 80 g
Excipients q.s.ad..... 100 g

Indications

Creatina Microsules improves muscular mass increase (weight) and thus, the animal's athletic performance. It increases energy, strength and vigor, especially in animal athlete which have delay fatigue. It improves post-effort recovery, achieving maximum strength and duration of the muscle contraction. It is recommended to obtain an optimal muscle profile in exhibit animals. The main effect of Creatine is the reversion (resynthesis) of burned energy (ADP) into available energy (ATP). After it is burned, the energy becomes ADP and it is renewed and transformed into ATP again (essential source of energy for the muscle contraction), this is why the animal athlete in training or in competition has less muscle fatigue, which will result in a longer-lasting energy reserve.

Target species

All animal species.

Route of administration

Oral; alone or mixed in the drinking water or in the feed. It is important this kind of product is administrated together with an additional carbohydrate (dextrose, glucose, etc.) for better absorption, apart from giving the animal an energetic plus. It is not recommended to stop the supplementation for more than 18 consecutive days. In case this happens, treatment should be resumed with the initial dose.

Dosage

* 0.24 g/Kg of animal weight.

Dose per 100 kg of l.w.: Initial dosage: 24 g every day. It may be divided into two or more applications during the day, for 14 consecutive days.
Maintenance dose: 8 g every day.

Precautions

It is not indicated in patients with renal insufficiency. Keep out of reach of children and pets. Never re-use empty containers; once used, packaging and surplus should be properly destroyed, either by burning or burying them in safe places away from lakes, streams, etc. Store in a cool dry place, away from sunlight.

Presentations

500 grs, 1 kg and 2 kgs.

Technical Director: Dr. Pablo Bringas
Made in Uruguay

Solo agent in U.A.E:
NATIONAL INDUSTRIAL COMPANY L.L.C
Tel: +971-2-6721000 – Fax: 971-2-6770011

Manufactured by:
LABORATORIOS MICROSULES URUGUAY S.A
Ruta 101 km 28, Camino Paso Escobar s/n
Pueblo Cap. Juan Antonio Artigas
Canelones, Uruguay. C.P.: 91.001
Tels: (00598-2) 2886761*, Fax: (00598-2) 2886760
www.laboratoriosmicrosules.com

